## Rasterelektronenmikroskopie

### Stichworte

Mikroskopie, Abbildung, Kryo, Zement, Hydratation, Mikrostructure

### Ansprechpartner

Günter Beuchle [guenter.beuchle@kit.edu](mailto:guenter.beuchle@kit.edu)

### Messprinzip

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## Scanning Electron Microscopy

### Keywords

Microscopy, Imaging, Cryo, Cement, Hydration, Micro Structure

### Contact

Günter Beuchle [guenter.beuchle@kit.edu](mailto:guenter.beuchle@kit.edu)


### Measuring principle

Scanning the sample surface point-by-point with an electron beam at vacuum conditions gives several imaging opportunities: every type of detector presents a different view of the scanned area.

Typically, pictures of signals from secondary electrons (SE) and backscattered electrons (BSE), respectively, are created. An energy dispersive X-ray detector is also available providing semi quantitative spot analyses of a given sample positions and even elemental mappings.

Microstructure imaging of cement water mixes can be achieved via a cryo preparation unit. Samples containing liquid water can be stabilized by rapid freezing. Water remains within the pore space without sublimation in the microscope vacuum. Thus changes of the microstructure of binder pastes as a function of contact – and hence hydration - time can be monitored.

The following figure shows an alite( $C_3S$ ) paste with a water to cement ratio of 0,5.

Samples were cryo prepared at three different hydration duration times. SE-images magnification increases from left to right. The densification of the microstructure is clearly visible.