
Röntgendiffraktometrie:

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X-ray diffractometry:

Quantitative phase analysis at room-temperature and *in-situ* measurements at temperatures up to 1200°C

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Method

X-ray diffractometry (XRD) is a method for investigation of various material properties, like the qualitative and quantitative phase composition or parameters of the crystal structures of individual phases. Apart from the chemical composition, these properties determine the behavior of materials. XRD utilizes the phenomenon of inelastic scattering of X-rays by periodic array of atoms within the crystal lattice. Thus, X-ray diffraction is limited essentially to crystalline samples with 3-dimensional ordered crystal structures.

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Application:

Phase identification and quantification of crystalline and nano-crystalline substances, determination of crystallinity, structure solution and refinement, non-ambient measurements, *in-situ* crystallization, recrystallization and reaction kinetics.

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Samples:

Powder, sample fragments, pastes, suspensions (slurries), thin layers etc. Amount of sample: 1mg-5g. Standard (e.g. Corundum) is added for quantitative analysis. This enables the determination of crystalline content and the calculation of amorphous content.

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Diffractometer:

Malvern –Panalytical Empyrean diffractometer. Cu tube, Bragg-Brentano HD optics, fixed slits, PIXcel 3D detector, soller slits (parallel collimators), divergence slit, sample changer for up to 45 samples.

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Measurement geometry:

Wide-angle X-ray scattering (WAXS) measurements in Bragg-Brentano $\theta/\theta\theta$ reflexion, transmission and small-angle x-ray scattering (SAXS).

Software: HighscorePlus, Easy SAXS, Topas 4.2, DiffracPlus, CrystalDiffract, Datenbanken: PDF2, COD19.

Software: HighscorePlus, Easy SAXS, Topas 4.2, DiffracPlus, CrystalDiffract, databases: PDF2, COD19.

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Sample holder: Backloading, frontloading, Si discs, sample carriers for transmission and for samples that are sensitive to air

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In-situ:

In addition to room temperature analytics, material properties can be investigated *in-situ* in a high-temperature chamber (HTK1200 Anton Paar). Maximum temperature is 1200°C and samples can be measured in air, N₂ and vacuum. Such measurements are used for monitoring of crystallization and re-crystallization processes, determination of thermal properties (expansion coefficient, density, structure parameter) and for the investigation of thermal-induced phase transformations.

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Room temperature measurements with standard provide comparable quantitative results. For the measurements in the high-temperature chamber, no standard is added in order to prohibit reactions between the standard substance and the sample.